

Evolution. Versatility. Simplicity.

Designed for Simple, Accurate Implantation and Reproducible Results

RESTORIS® MetaFix™ Hip Stem



Used with the RIO® Robotic
Arm Interactive Orthopedic
System to enable a new
level of precision


MAKO
SURGICAL CORP.®

RESTORIS® MetaFix™ Hip Stem

A collarless, cementless, fully coated hydroxyapatite implant available in a broad range of sizes, allowing surgeons to choose the right solution for individual patient's needs.

RESTORIS® MetaFix™ is an optimal solution for a broad range of hip defects, available with two different neck offset options and sizes that match a wide variety of anatomy.



Clinically Proven Design¹

Dual-taper trapezoidal design with vertical and horizontal grooves which may reduce rotational and axial displacement and may promote stability.



Versatility

Two neck offset options allow for the natural restoration of patient's biomechanics and provide an optimal solution for a broad range of hip defects.

Bone-conserving Design

Low profile shoulder is bone conserving in the greater trochanteric region and allows for ease of use in smaller incision surgery.



Hydroxyapatite Coating

Biocompatible hydroxyapatite coating thickness of 155µm is designed to provide initial press-fit and may help to promote fixation².



Simplicity

Designed to be ergonomic and provide procedure efficiency, the RESTORIS® MetaFix™ instrumentation is suitable for all surgical approaches.



1. Hallan G, Lie SA, Furnes O, Engesaeter LB, Vollset SE, Havelin LI. Medium and long term performance of 11516 uncemented femoral primary stems from the Norwegian arthroplasty register. J Bone Joint Surg [Br] 2007 Dec;89(12):1574-80.

2. Hardy D, Frayssinet P. Hydroxyapatite-coated femoral arthroplasties: A long term study through 29 Corail prostheses explanted during a ten-year survey. Surgical Technology International X. 2003.



Restoring Quality of Life Through Innovation®

MAKO Surgical Corp.
2555 Davie Road | Fort Lauderdale, FL 33317 | 866.647.6256 | www.makosurgical.com

Manufactured by Corin Group PLC.
The Corinium Centre, Cirencester, GL71YJ, United Kingdom
Phone: 011 44 1285659866

Distributed by MAKO Surgical Corp.

All third party trademarks used herein are property of their respective owners.

© MAKO Surgical Corp. 204937 r00 05/11